Outside

MADE WITHOUT
BPA/BPS

Hot & Cold Breast Therapy Packs

Cold: Reduce pain and swelling due to engorgement.
Cold: Relieve breasts with plugged ducts or mastitis.

Hot: Encourage milk flow.

Contents: 2 reusable packs with soft covers

INTENDED USE:

Cold Therapy can help ease the symptoms of engorgement. Engorgement occurs when the breasts become overfull and they become swollen and painful. It can happen at the beginning of breastfeeding when your mature milk comes in or at the end of breastfeeding when you are weaning. Apply cold packs between feeding to help ease swelling and pain.

Cold Therapy can help ease plugged ducts and mastitis. Plugged ducts happen when there is a buildup of milk solids in the ducts of your breasts and it prevents milk from flowing through. Mastitis is when you get an infection in combination with a plugged duct. If you experience flu-like symptoms or have a fever, be sure to reach out to your healthcare provider. Applying cold therapy before feeding or pumping along with lymphatic breast massage, can help break up the blockage.

Hot Therapy can help encourage milk let-down and make pumping more effective. Apply warm packs before your pumping session or wrap the warm packs around standard pump flanges to reduce the amount of time spent pumping.

DIRECTIONS:

Lansinoh® Hot & Cold Breast Therapy Packs may be applied for up to 20 minutes, or as doctor recommended.

For cold therapy:

Freeze for at least 4 hours.

For heat therapy:

Place one pack inside the microwave lying flat for time indicated in the chart. Each pack should be heated separately. After testing pack temperature, you may


Suggested heating time for (1) Hot & Cold Breast Therapy Pack	
700 Watts	16 Seconds
1,000 Watts	13 Seconds
1,250 Watts	11 Seconds

add more time in 5 second increments as needed until desired temperature is reached. Always test pack temperature prior to application.

IMPORTANT: Do not heat for more than time indicated below during initial heating and no more than 5 second increments during subsequent heating. DO NOT HEAT MORE THAN 25 SECONDS TOTAL.


Use without pump:

Cool or heat as directed and place inside the soft cover prior to placing on the breast. The nipple is located at the center of the pack. To hold them in place, packs can be placed in bra. Do not use longer than 20 minutes in a single session.

Use with pump:

Compatible with manual and standard electric pumps. Heat as directed and put the pack within the soft cover around the breast pump's flange. Lock the soft cover's ends with the press button. The pack stays around the flange during pumping. Do not use longer than 20 minutes in a single session.

Cleaning:

After each use, clean the pack with mild soap and water or 70% rubbing alcohol. We do not recommend the use of ethanol, phenol, or chlorine-based products.

Store in freezer or airtight container, out of direct light.

CAUTION: Always test pack temperature prior to application.

Prolonged application may cause injury, including burns. Do not use on infants. Do not use if you have poor circulation. Do not use heat therapy on stitches or an open wound. Do not fall asleep using pack. Use only under adult supervision. For extreme, prolonged, or chronic pain, contact a physician. Use only as directed. Heating beyond suggested time may damage or destroy the pack. Do not puncture and/or ingest contents. For external use only. Ensure pack is clean before use. Do not place soft covers in microwave.

Distributed by: / Distribuido por:
Lansinoh Laboratories, Inc., Alexandria, VA 22314, USA | +1 800-292-4794 | Lansinoh.com
Made in China / Hecho en China

Lansinoh® is a registered trademark of Lansinoh Laboratories, Inc. in the United States. / Lansinoh® es una marca registrada de Lansinoh Laboratories, Inc en los Estados Unidos. © 2024 Lansinoh Laboratories, Inc.

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All rights reserved / Reservados todos los derechos.